**WAGNERSIL
Gebrauchsinformation**

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WAGNERSIL
Instructions for use

WAGNERSIL is an addition-curing 1: 1 duplicating silicone (RTV 2K silicone rubber) with different final hardnesses and processing and curing times. It is low-viscosity with excellent mechanical properties. WAGNERSIL is used to produce silicone impressions that are extremely accurate in terms of dimensions and details. In its final consistency, WAGNERSIL has a very high tear strength and is very flexible. The model can therefore be easily removed from the mold, especially with undercuts. WAGNERSIL does not require any wetting agents. As a result of the addition cross-linking with a platinum catalyst in the A component, the silicone hardens without shrinking and odorless. WAGNERSIL is also very suitable for the flask-free technique. WAGNERSIL can be mixed simply and easily both in a weight ratio and in a volume ratio of 1: 1. WAGNERSIL can be poured several times with almost any casting compound. It is temperature resistant up to 200 ° C, briefly up to 250 ° C. WAGNERSIL is ideally suited for all impressions / molds in the fields of industry, craft and hobby. All WAGNERSIL silicones are physiologically harmless (can be used on the skin). The silicone types WAGNERSIL 26 LE, 28 LE and 35 LE are food-safe duplicating silicones. The food-safe silicones meet the requirements of the BfR and FDA.

Processing information:
Wagnersil (all variants)

Place the model cleaned with alcohol or white spirit in a suitable mold or cuvette. Then mix WAGNERSIL components A and B at room temperature in a ratio of 1: 1 for approx. 1 minute until a streak-free, uniform shade is obtained. Then pour the silicone into the same place with a thin stream until the model is covered all around. Processing time see table. Release agent (Wagnersil TMV) is required when silicone is poured onto silicone or when using special model materials. At a room temperature of 23 ° C, the model can be removed after the curing time (see table for demolding time). The curing time can be shortened by adding heat (e.g. at 50 ° C). Only pour out the silicone mold after a further 20 minutes so that any deformations that may have arisen from demolding can be reset.

Note: If the material or room temperature is above 23 ° C, the processing time is significantly reduced.

Attention: Some model materials contain "platinum inhibitors". These prevent the silicone from hardening on the contact surface with the model.

Recommendation: clean the model well and do a test with a small amount before pouring, if necessary isolate the model with Wagnersil TMV.

WAGNERSIL food safe

Variants: Wagnersil 26 LE, WAGNERSIL 28 LE, WAGNERSIL 35 LE

(Fulfill the BfR and FDA requirements for food-safe duplicating silicone)

When used with food, the silicone mold must be degassed or tempered. Place the mold in the oven for approx. 4 hours at 200 ° C

Technische Daten / Technical data:

Produkt Daten / Typ	9N	17N	19NO	20NF	21NO	22N	22S	22NF	25NF	25L	26LE	28LE	32N	35LE
Endhärte, Hardness (Shore A)	9	17	19	20	21	22	22	22	25	25	26	28	32	35
Verarbeitungszeit, Pot life (min.)	6	6	5	10	6	6	3	6	6	30	30	15	5	6
Entformzeit, demoldable after (min.)	30	30	15	90	30	30	12	30	30	360	360	60	30	60
Viskosität, Viscosity (mPas)	2000	2400	9000	4500	2500	2500	2500	5000	5500	5000	5200	6000	9500	6000
Dichte, Density (g/cm³)	1,05	1,05	0,96	1,12	1,05	1,05	1,05	1,13	1,14	1,05	1,10	1,05	1,19	1,10
Reißdehnung, Elongation at break (%)	260	320	180	300	410	410	410	330	480	500	430	500	380	450
Weiterreißfestigkeit Tear strength (N/mm)	2,4	3,8	2,5	3,2	5,0	5,0	5,0	3,6	7,5	25,0	11,0	25,0	8,3	14,0
Reißfestigkeit, Tensile strength (N/mm²)	1,5	3,2	2,5	2,9	3,1	3,1	3,1	3,9	4,0	5,0	4,5	5,0	4,3	5,0
Farbe, Color	Blau blue	Blau blue	Klar clear	Grün green	Transluzent translucent	Rosa pink	Rosa pink	Rosa pink	Rosa pink	Flieder lilac	Lavendel lavender	Weiß white	Blau blue	Transluzent translucent